ZEISS

West Germany

Carl Zeiss

D-7082 Oberkochen

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
Tab. 1: Strahlenteilung an beiden Tuben
Beam splitting in both tubes

U	B	O	A1	A2
1	100%	1,2,3	0%	0%
2	20%	1	0%	80%
2	20%	2	40%	40%
2	20%	3	80%	0%
3	0%	1	0%	100%
3	0%	2	50%	50%
3	0%	3	100%	0%

Installation of the Axiophot camera system on the Axioplan microscope

The Axiophot (45 19 37-9902) and Axiophot Pol (45 19 36-9902) camera systems can be mounted on the Axioplan microscope (from April 1990) as follows:

- Remove tube panel and binocular tube from the Axioplan stand. Screw out both pins (4) using a screwdriver.
- Screw out retaining screw (1) until it is no longer visible in the dovetail mount.
- Lift camera system by the sides with both hands. Insert camera system in the tube mount via the dovetail and align it parallel to the stand edge.
- Move camera system until the thread (3) in the Axioplan stand becomes visible.
- Insert screw and washer in the camera system and the thread of the Axioplan stand, but do not tighten it fully.
- Tighten tube retaining screw (1).
- Tighten back screw (2).
- Plug the connecting cable into the camera system and into the camera control panel and screw the connectors tightly in position.

Caution: Plug in or remove cable only if instruments are disconnected from the power line.

Components and operation of the module are described in:
 G 42-250/II Axiophot Transmitted Light and Incident-Light Fluorescence
 G 42-350/II Axiophot Reflected light
 G 42-351 Axiophot Pol Transmitted and Reflected light

Installation of the bin. photo/TV tube with 2 ports on the Axioplan microscope

This module is installed in the same way as the Axiophot camera system.

Photo and TV tubes can be attached to the two ports, as shown in the illustration.

The two pushrods split the beam as follows (see also Tab.1):

<u>lower pushrod U</u> pushed in (U1):	Observation only (100% of the light to the tube).
pulled out to 1st stop: (U2)	20% of the light for observation, 80% to upper beam splitter.
pulled out all the way: (U3)	no observation through the tube no light to the eyepiece, 100% to upper beam splitter.
<u>upper pushrod O</u> pushed in (O1):	100% of the light to port A2.
pulled out to 1st stop: (O2)	50% to port A1 50% to port A2
pulled out all the way (O3):	100% of the light to port A1.

The light intensity at rear port A2 can be attenuated using one of the three neutral density filters supplied (0.6, 0.3 and 0.1) to ensure that it is in compliance with the sensitivity of the different TV cameras. (Remove mount from A2 completely, replace plane glass with one of the filters and insert mount again.)

2